

# Neater Eater Robot



## Overview

The Neater Eater Robotic enables you to feed yourself. It is simple to use; controlled by either a Tablet or plug-in switches.

Tablet menus and help guides demonstrate many options such as setting your mouth position, automatic spoon-wiping and plate-turning. There are fork and snack-holder modes and a routine to scrape the plate at the end of your meal.

- (A) On/Off switch
- (B) Removable touchscreen Tablet with onscreen help and video guides
- (C) Robot arm has wide range of movement to reach your mouth
- (D) Choice of plug-on cutlery
- (E) Choice of plates that rotate


*PTO for easy setup!*

1 Turn ROBOT on

2 Turn TABLET on and wait for onscreen guide

😊 Enjoy your meal!

*(Remember to charge the Robot and Tablet every night)*



Removable touchscreen


Tablet for easy access.

Replace here between meals to recharge its battery wirelessly.

Red Light = charging  
(flashing = battery low)  
Green Light = fully charged

Plug power in here

overnight to charge batteries  
(can be plugged in when eating).



Orange Light = Robot battery charging  
Red Light = power to Tablet charger  
Green Light = external power plugged in  
Blue Light = Tablet Bluetooth connected

### Problem?

Ensure robot and tablet are charged.  
Hold Tablet power and the volume up buttons until it vibrates.  
Hold Robot power button until light goes out.  
Press Robot power button again—it lights up.  
Wait two minutes for Neater App to start.



Press this question mark icon on screen for help and hold for videos.

The Tablet will guide you through set-up and use.  
PLEASE READ ALL SAFETY INSTRUCTIONS.

For further information, video hints, and tips go to [www.neater.co.uk](http://www.neater.co.uk)