Thera Standard/150 Manual height adjustment



Thera Gas/150
Gas spring height adjustment



Thera Double/130 Electrical hight adjustment and electrical leg spreading



Thera El/150 Electrical height adjustment



Thera Extra/240 Max. capacity up to 240kg

# Thera

Walking Tables: Safe and stable support aid

# Primary application

Home care Institutional Hospital

The Thera Walking Table programme consists of a series of safe and stable ambulatory devices.

They are designed to support the weight of the user when gait training in a physio or rehab situation, or in a post operative situation by encouraging early patient mobility.

These walking tables can be used independently or with assisted care. They are easily height adjusted to accommodate users of different heights.

The range consists of 5 models.

### Key features:

Rear braked wheels Width adjustable arm padding and handgrips Manual, gas or electrical height adjustment Stable non-tip base Maximum weight capacity of 130, 150 or 240 kg

#### Accessories:

Dual hand brakes Side support Foot plate Dripstand holder





Dual hand brakes accessory.



Multi adjustable handles to suit varying needs. Width adjustable armpads.




Thera Extra: Double columns for increased bariatric stability.



Thera Gas: Easy height adjustability using the gas spring release lever.



Thera Double has electrical leg spreading for extra space around the client.



## **TECHNICAL SPECIFICATION**

#### THERA STANDARD

Max. load150 kgWidth outer780 mmHeight adjustement900-1320 mmWeight of Thera13 kg

**THERA GAS** 

Max. load150 kgWidth outer780 mmHeight adjustement900-1290 mmWeight of Thera15 kgHeight of Thera910-1260 mm

THERA ELECTRIC

Max. load150 kgWidth outer780 mmHeight adjustement900-1300 mmWeight of Thera21kg

THERA DOUBLE

Max. load130 kgWidth outer640-850 mmHeight adjustement870-1280 mmWeight of Thera31 kg

THERA EXTRA

Max. load240 kgWidth outer880 mmHeight adjustement920-1320 mm

Weight of Thera 35 kg

Consult the Technical manual for further specifications. All figures and measurements are for guidance only and may differ in actual situations. Ergolet reserves the right to changes without notice.

Product photos may include accessories. Consult your dealer.



Ergolet
Taarnborgvej 12C
4220 Korsoer
Denmark
tel. +45 70 27 37 20
e-mail: ergolet@ergolet.dk

